

3pp[®] Radial Hinged Ulnar Deviation[™] Splint




The Radial Hinged Ulnar Deviation Splint aligns each finger individually to **improve the alignment of strongly deviated fingers**. The radial hinge allows free metacarpophalangeal (MP) flexion and extension for improved hand function.

Features and Benefits

The Kydex plastic shell is lined in soft felt and can be adjusted with a heat gun for a custom fit.

Finger strap allows each finger to be adjusted individually for proper alignment.



The strong radial hinge aligns with the axis of motion at the 2nd MP joint for motion.



Open the finger straps and loosen the palmar stay for easy application.

The hypothenar extension and the palmar stay help secure the brace on the hand.



Indications

Conditions resulting in ulnar or radial deviation at the MP joints:

- Rheumatoid Arthritis (adult and Juvenile Idiopathic Arthritis)
- Lupus
- Neurological Conditions

Suggested Billing Code

HCPCS L3929 - Use of this code does not guarantee reimbursement.

Sizing Information

Lay hand on a flat surface and measure across the MP joints of digits 2 through 5.

Size	Fits
Small	2 3/4" to 3"
Medium	3" to 3 3/8"
Large	3 3/8" to 3 5/8"

