





Oval-8[®] CLEAR Finger Splints provide the same easy-to-wear, comfortable support and features as the classic beige Oval-8 splints. Oval-8 CLEAR splints blend in with all skin tones and can be turned on the finger to keep **mallet** or **jammed fingers** straight; limit motion to stop **trigger finger** and **trigger thumb**; control **hypermobility** and protect and rest **crooked** fingers.

Each splint helps manage 6 or more conditions





Trigger Finger





A Simple Turn of the Splint Changes How an Oval-8 Fits

Put the end with the plus + sign on first for a looser fit or turn it around for a tighter fit.



Find the 🕇 sign on the top of the splint





Customize Oval-8 Splints in under 20 seconds



Oval-8 splints can be adjusted using a narrow nozzle heat gun to alter fit, length and the angle the splint holds the finger. CLEAR splints require less heating time than the Beige splints and a more delicate touch when making adjustments.



Oval-8 CLEAR Finger Splint FAQs

Are the Oval-8 CLEAR Splints the same as the BEIGE Splints?

Both the Classic Beige splints and the NEW CLEAR splints are made in the same molds but they are made from different materials.

- · Beige splints are 100% Polypropylene and CLEAR splints are made from Copolyester
- There may be minor variations in sizing but adjustments can be made to assure a proper fit The CLEAR splints are more resistant to bending and it is best to avoid repeatedly stressing the splints

Can Oval-8 Splints be adjusted?

Both the Beige and Oval-8 CLEAR splints can be adjusted with a heat gun with a narrow nozzle
CLEAR splints require less heating time (2-3 seconds) than the Beige (3-4 seconds) and require less pressure to alter the splint

How long will Oval-8 Splints last?

How the splints are worn and used will determine how long the splints last. The Beige splints have been proven to last for several years of daily use. It is not yet known how long the CLEAR splints will last.

Do they contain latex?

No natural rubber latex was used in the product or packaging.

Ways to Purchase Oval-8 CLEAR Finger Splints

Sizing Set



With one of each of the 14 Oval-8 sizes, the Sizing Set is the most accurate way to find the correct size for patients in seconds.



Graduated Sets with 3 sizes and Combo Sets with 5 sizes provide a mix of sizes for patients whose finger may change over time.

Single Size Packages



Oval-8 splints are available in single size packages of 1 or 5 splints to issue to patients or refill an Oval-8 Kit. *Kit is not available with CLEAR splints.*

3point: Health Care Providers Customer Service

888-378-7763 Mon-Fri 8:30am - 4:30pm EST service@3pointproducts.com

ohmyarthritis

Patients & Consumers Online Ordering

www.shopoma.com

Convenient, cost effective ordering