

3pp® Ez FIT ThumSpica Splint



- Carefully cushions with a tri-layer memory foam to fully support the thumb and forearm
- Full length aluminum stay can be custom contoured to position the thumb and wrist
- Fits along the radial side of the forearm, leaving fingers and the ulnar side of the hand free for function
- Flat buckles prevent straps from twisting and straps are angled to contour to the arm for a secure fit
- One size, fits left or right hand for economic stocking

Indications

- de Quervain's Tendinosis
- Skier's / Gamekeeper's Thumb
- Overuse / Repetitive Stress
- Arthritis

Features and Benefits



Sizing Information



Sizing Information
Measure around your wrist and mid-forearm. Fits both right/left hand.

One size fits most
Wrist:
4.5" to 8.5" (11.4 - 21.6cm)
Mid-forearm:
Fits up to 10.5" (26.7cm)

Healthcare Professionals

888-378-7763 M-F 8:30 - 5:00 EST
service@3pointproducts.com

Order Number
One Size fits both Right & Left

One Size	P3014-U
----------	---------

Patients & Consumers

Online Ordering
www.ohmyarthritis.com



[ohmyarthritis.com](http://www.ohmyarthritis.com)