

3pp[®] Carpal Lift[™] Brace



The 3pp Carpal Lift employs a unique strap & pad system that applies counterforce support to reduce wrist pain due to ulnar or mid-carpal instability. A padded pocket connects to adjustable straps that thread through the splint to lift and correctly position the wrist. The lifting motion lessens the clicking and popping sensation felt when gripping or turning the palm and reduces stress on the carpal ligaments and TFCC.

Features and Benefits



Adjustable straps thread through the brace to lift and position the wrist to reduce pain from the popping and clicking felt when gripping or turning the palm



The palmar pocket features a cushioned pad and adjustable straps



Pad applies a corrective lift of the ulnar and mid-carpals and reduces stress on the TFCC



Press down on the palm with the wrist in neutral for optimal lift using the adjustable straps



If needed, small slits in the thumb hole will help position the splint for a custom fit

Indications

- Ulnar sided wrist pain
- Ulnar carpal instability
- Mid carpal instability
- TFCC injury
- Guyon's Canal syndrome

Suggested Billing Code

HCPCS L3908 - Use of this code does not guarantee reimbursement.

Sizing Information

Choose right or left hand
Measure around the hand at the big knuckles (MP joints)

Sml/Med 6" to 7¹/₄" (15.2 cm - 18.4 cm)
Med/Lrg 7¹/₄" to 9" (18.4 cm - 22.8 cm)

Hint: If you are at the end of a size range

- Go up in size if your hand is muscular
- Go down in size if your hand is thin

