



3pp® Buddy Loops® & Toe Loops®

3pp® Buddy Loops®



Non-slip, foam-lined 3pp Buddy Loops are the most effective way to secure fingers to protect or assist motion for **jammed or sprained fingers**.

1/2" (1.27 cm), 3/4" (1.9 cm) & 1" (2.5 cm)
 All sizes are 5" (12.7 cm) long




3pp® Toe Loops®



The adjustable 3pp Toe Loops provide cushioned relief for **Hammertoes, Claw Toes, Overlapping Toes, and Fractures**.

Narrow measures 7/8" wide and 1 5/8" long. Wide measures 1 1/4" wide and is 1 1/8" long.



FAQs

3pp Buddy Loops

What is the difference between the 3 sizes?

- 1/2" is best for smaller hands and light activities
- 3/4" is a good all around size
- 1" is best for large hands and extra support

3pp Toe Loops

How do I choose Narrow or Wide?

- Narrow generally fits the four small toes
- Wide is designed to wrap the big toe and the first small toe for support and following surgery

To Order from 3-Point Products

Health Care Providers

888-378-7763

Mon-Fri 8:30am - 4:30pm EST

service@3pointproducts.com

Online Ordering

Patients & Consumers

www.ohmyarthritis.com



Scan Here