

# Fix Comfort Thumb Brace

Manutec® by Orliman®



- Tri-layer of foam and velour material cushions the brace
- Aluminum interior frame contours to the hand and thumb
- Soft edges eliminate pressure around the wrist and hand, and perforations keep skin cool and dry
- Universal fits right & left in Sm/Med and Med/Lrg

## Indications

- de Quervain's Tendinosis
- Gamekeeper's / Skier's Thumb
- Ligament Injuries
- CMC Thumb Arthritis
- Sprains, Strains

## Features and Benefits




Lightweight aluminum contours for custom fit



Fits right/left hand



Adjustable thumb strap for optimal comfort



3-Layer breathable material keeps skin cool and dry



## Sizing Information



Measure around your wrist to find the best size for you

Small/Med  
5" to 7" (14 - 18cm)

Med/Large  
7" to 8¾" (18 - 22cm)

## To Order from 3-Point Products

Health Care Providers

888-378-7763 Mon-Fri 8:30am - 5:00pm EST

customerservice@3pointproducts.com

## Online Ordering

Patients & Consumers

[www.ohmyarthritis.com](http://www.ohmyarthritis.com)



[ohmyarthritis.com](http://ohmyarthritis.com)

## Order Numbers

Sm/Med	Med/Lrg
P2020-23	P2020-34



flip the page